



SPECIFICATION

Measuring Modes: body temperature
milk temperature

Non-touch: Avoid contact with human skin to avoid cross infection

Measuring Distance: 1-5cm

Allowance Difference:

42.1°C-42.9°C, $\pm 0.3^{\circ}\text{C}$ (body temperature)

35.0°C-42.0°C, $\pm 0.2^{\circ}\text{C}$ (body temperature)

32.0°C-34.9°C, $\pm 0.3^{\circ}\text{C}$ (body temperature)

Measuring Range:

32.0°C-42.9°C (body temperature)

10.0°C-60.0°C (milk temperature)

Display: LCD

Low Battery Reminder: $<2.6\pm 0.1\text{V}$

Display Resolution: 0.1°C

Unit: Celsius degree, Fahrenheit degree

Working Circumstance: 10°C-40°C

Storage Circumstance: -20°C - $+55^{\circ}\text{C}$

Power Saving Mode: power off automatically if no operate after 1 min

Memory Record: 25 records

Battery: 2 x AA

Parating Life: 2 years

Product Size: 161 x 79 x 46mm

GUANG DONG JU FENG PRECISION TECHNOLOGY CO.,LTD

Factory address:

Fuxing Road, No.2 Tianwei Industrial Park,

Changping Town, Dongguan City, Guangdong

Province, China,

Tel: +86-769-85328518 86997566

Fax: +86-769-85414142

http://www.dgjufeng.com

DIGITAL INFRARED THERMOMETER

Safety • Fast • Precised • Non-touch • Multi functional



JF-TWT201

- **BODY TEMPERATURE MEASURING**
- **MILK TEMPERATURE MEASURING**
- **BABY BATH WATER TEMPERATURE MEASURING**
- **BABY ROOM TEMPERATURE MEASURING**

DIGITAL INFRARED THERMOMETER

1. Body Temperature Mode:

measuring range between 32-42.9°C (89.6°F-109.22°F)

BODY TEMPERATURE



2. Room Temperature Mode: 0-40°C (32°F-104°F)

MILK TEMPERATURE



3. M.K Temperature Mode:

Milk temperature measurements range from 10 to 60 degrees (50-140°F) and can measure the temperature of any object with a temperature without touching it

Baby Bath Water:

When measuring the temperature, please measure it from directly above the water. It is recommended that if the bath water is left for too long, it is better to shake the water well first to ensure the accuracy of the average temperature.

It's suggested that your baby's bath water be between 37 and 43 degrees



Milk Temperature Measuring

When measuring the temperature, shake well and then measure the water level of the bottle.

It's recommend to the baby to eat milk temperature control in 37 degrees -38 degrees, so that the baby's stomach function development better, more easily absorb and digest nutrients

